


Design Matters Competitive Comparison

Not all matrices are created equal

In benchtop testing, Integra DRT, PriMatrix and PolyNovo NovoSorb® Biodegradable Temporizing Matrix (BTM) were evaluated to compare physical and biological attributes.¹

Cell adhesion, migration, and proliferation are required for robust and efficient wound healing.² Physical properties of a matrix, such as hydrophilicity, conformability, and porosity are important factors to aid in the healing process.³



Matrix Comparison

Integra® DRT

Collagen and Chondroitin-6-Sulfate
with Silicone Layer



Silicone:

Protects the Wounds

- Provides immediate coverage, acting as temporary epidermal layer⁴
- Maintains a moist wound environment⁵
- Typically removed between 14-21 days

Collagen and Chondroitin-6-Sulfate: Promotes Dermal Regeneration

- Bioengineered acellular matrix minimizes inflammatory response^{6,7}
- 3D pore structure optimized for cellular and vascular ingrowth^{8,9,10}
- Grows with the patient and helps restore function and joint mobility^{6,11}

PriMatrix®

Bovine-derived
Dermal Repair Scaffold




Acellular Fetal Bovine Dermis

- Collagen I and III are associated with wound healing and developing tissues^{12,13,14}
- Robust structure for enhanced handling and tensile strength^{13,15,16}
- Supports revascularization and re-epithelialization processes^{17,18,19}

NovoSorb® BTM

Biodegradable Polyurethane
with Sealing Membrane



Polyurethane Sealing Membrane

- Temporary non-biodegradable layer limiting moisture loss
- Serves as a barrier to outside bacteria and closing wound^{20,21}
- Typically removed between 3-5 weeks (i.e., 21-35 days)

Biodegradable Polyurethane Matrix

- Designed as a temporizing matrix
- Open cell structure with large chambers and pores
- Rigid mechanical properties


Adherence and Proliferation^{4,21,22,23}

When compared to NovoSorb® BTM, cells seeded in Integra® and PriMatrix® demonstrated:


- ✓ **Over 200% faster** fibroblast adhesion
- ✓ **Over 2X** sustained fibroblast adherence levels
- ✓ **Over 2X** (Integra) and **over 1.5X** (PriMatrix) fibroblast proliferation
- ✓ **Over 20X** (Integra) and **over 10X** (PriMatrix) adherent keratinocytes

Note: Preclinical results are not necessarily indicative of clinical outcomes.

Adherent Fibroblasts



Adherent Keratinocytes



Physical Comparison

Not all matrices are created equal


Porosity

The overall porosity and size distribution of pores may impact how cells migrate into a matrix to support a variety of essential cellular activities.⁵

Integra's bioengineered matrix and PriMatrix's dermal scaffold align to native human architecture with small, tight pore distribution allowing for cell migration, adherence, and proliferation.¹⁴

Integra® DRT


Collagen-C6S
with Silicone Layer



Engineered design offers optimized pore size and structure, allowing cell migration, repopulation, and remodeling¹⁴

PriMatrix®


Bovine-derived
Dermal Repair Scaffold



Biocompatible porous matrix, which provides the foundation for rapid cellular repopulation and tissue revascularization to support the wound healing process^{13,18,19}

NovoSorb® BTM

Biodegradable
Polyurethane




Wide range of larger pore sizes less ideal for cell migration, repopulation, and remodeling¹

Absorption


Hydrophilic materials have an environment conducive for autologous blood, which contain essential cells and growth factors for infiltration, attachment, and repopulation.

Hydrophobic materials create a challenging environment for blood to absorb and spread throughout the matrix.^{1,4}


- Blood absorbed into Integra and PriMatrix exemplifies resilient hydrophilic properties.^{1,22}
- BTM repelled blood, exhibiting hydrophobic properties.^{1,22}



Bioengineered 3D matrix with an optimized design to help promote cellular ingrowth



Blood and cells are absorbed, providing a foundation for cellular repopulation and tissue revascularization^{4,5}



Exhibits prolonged hydrophobic properties with blood beading on top of the matrix, resisting absorption

Conformability

High conformability allows matrices to contour to complex wound geometries for intimate contact, permitting cells to migrate, repopulate, and remodel. Both Integra and PriMatrix are substantially more conformable than BTM.^{1,22}



Integra is approximately 9x more comfortable than BTM



PriMatrix is approximately 2X more conformable than BTM



BTM provides minimal conformability

Note: Preclinical results are not necessarily indicative of clinical outcomes.

Brief Summary

Consult Package Insert for Full Information

Integra® Dermal Regeneration Template

Indications

Integra template is indicated for the postexcisional treatment of life-threatening full-thickness or deep partial-thickness thermal injuries where sufficient autograft is not available at the time of excision or not desirable due to the physiological condition of the patient. Integra template is also indicated for the repair of scar contractures when other therapies have failed or when donor sites for repair are not sufficient or desirable due to the physiological condition of the patient.

Warnings

Hemostasis must be achieved prior to applying Integra template. Inadequate control of bleeding will interfere with the incorporation of Integra template.

PriMatrix® Dermal Repair Scaffold

Indications

PriMatrix is intended for the management of wounds that include partial and full thickness wounds, pressure, diabetic, and venous ulcers, Second-degree burns, surgical wounds—donor sites/grafts, post-Moh's surgery, post-laser surgery, podiatric, wound dehiscence, Trauma wounds—abrasions, lacerations and skin tears, tunneled/undermined wounds, Draining wounds.

References:

1. Data on file at Integra LifeSciences Corporation: REGEN-21-004. 2. Wilkinson, H. N., & Hardman, M. J. (2021). Wound healing: cellular mechanisms and pathological outcomes. *Open Biology*, 10(9), 200223. 3. Anselme K, Ploux L, Ponche A. Cell/Material Interfaces: Influence of Surface Chemistry and Surface Topography on Cell Adhesion. *J Adhes Sci Technol*.2010;24(5):831-852. 4. Burke JF, Yannas IV, Quinby WC, et al: Successful use of a physiologically acceptable artificial skin in the treatment of extensive burn injury. *Ann Surg*. 1981 (Oct);194(4):413-428. 5. Bryers, J. D., Giachelli, C. M., & Ratner, B. D. (2012). Engineering Biomaterials to Integrate and Heal: The Biocompatibility Paradigm Shifts. *Biotechnology and Bioengineering*, 109(8), 1898-1911. 6. Burke JF: Observations on the development and clinical use of artificial skin: An attempt to employ regeneration rather than scar formation in wound healing. *Jpn J Surg*. 1987;17(6) 431-438. 7. Orgill DP, Butler CE, Regan JF. Behaviour of collagen-GAG matrices as dermal replacement in rodent and porcine models. *Wounds: A Comp of Clin. Research and Practice*. 1996;8(5):151-157. 8. Boo S, Dagnino L. Integrins as Modulators of Transforming Growth Factor Beta Signaling in Dermal Fibroblasts During Skin Regeneration After Injury. *Adv Wound Care (New Rochelle)*. 2013 Jun;2(5):238-246. 9. Boekema BK, Vlieg M, Olde Damink L, Middelkoop E, Eummelen L, Bühren AV, Ulrich MM. Effect of pore size and cross-linking of a novel collagen-elastin dermal substitute on wound healing. *J Mater Sci Mater Med*. 2014 eb;25(2):423-433. 10. Später T, Menger M, Nickels R, Menger M, Laschke M. Macrophages promote network formation and maturation of transplanted adipose tissue-derived microvascular fragments. *J Tissue Eng*. 2020 Feb;11:1-14. 11. Sheridan RL, Hegarty M, Tompkins RG, Burke JF: Artificial skin in massive burns—results to ten years. *Eur J Plastic Surgery*. 1994;17:91-93. 12. Berger A, Tanzella U, Machens HG, Liebau J. The use of Integra in primary burn wounds and unstable secondary scars. *Chirurg*. 2000; 71(5):558-563. 13. Cornwell KG, Landsman A, James KS. Extracellular matrix biomaterials for soft tissue repair. *Clin Podiatr Med Surg*. 2009 Oct;26(4):507-23. 14. Mays PK, Bishop JE, Laurent GJ. Age-related changes in the proportion of types I and III collagen. *Mech Ageing Dev*. 1988 Nov 30;45(3):203-12. 15. Ramshaw JA. Distribution of type III collagen in bovine skin of various ages. *Connect Tissue Res*. 1986;14(4):307-14. 16. Data on File with Integra LifeSciences. 17. Neill et al.Utilizing Biologic Assimilation of Bovine Fetal Collagen in Staged Skin Grafting, *Annals of Plastic Surgery*, 2012. 18. Rennert R, et al. Cellular Response to a Novel Fetal Accellular Collagen Matrix: Implications for Tissue Regeneration. *International Journal of Biomaterials*. 2013. 19. Lullove, E. Accellular Fetal Bovine Dermal Matrix in the Treatment of Nonhealing Wounds in Patients with Complex Comorbidities. *Journal of the American Podiatric Medical Association*. 2012; 102 (3): 233-239. 20. Greenwood JE, Dearman BL. Comparison of a sealed, polymer foam biodegradable temporising matrix against Integra(R) dermal regeneration template in a porcine wound model. *J Burn Care Res*. 2012; 33:163-73. 21. Wagstaff MJD, Schmitt B, Caplash Y, Greenwood JE. Free flap donor site reconstruction: A prospective case series using an optimized polyurethane temporising matrix. *Eplasty*. 2015; 15:231-48. 22. Stefanelli VL, Mintz B, Gandhi A, Smith J. Design matters: a comparison of natural versus synthetic skin substitutes across benchtop and porcine wound healing metrics: an experimental study. *Health Sci Rep*. 2023;6:e1462. doi:10.1002/hsr2.1462. 23. Proceedings of the National Academy of Sciences USA. Vol. 86, Feb. 1989, pp. 933-937.

A surgeon must always refer to the package insert, product label and/or instructions for use before using any Integra LifeScience product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Integra LifeSciences representative if you have questions about the availability of products in your area.

Integra, PriMatrix, and the Integra logo are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries. All other trademarks are trademarks of their respective owners or holders.

Ordering Information

Not all size options and configurations are shown

Product Size		Solid	Meshed
IDRT	2x2in	32021	MDRT2021
	4x5in	34051	MDRT4051
	4x10in	34101	MDRT4101
	8x10in	38101	MDRT8101
Bilayer	2x2in	BMW2021	MWM2021
	4x5in	BMW4051	MWM4051
	4x10in	BMW4101	MWM4101
	8x10in	BMW8101	MWM8101
IWM	2x2in	52021	
	4x5in	54051	
	4x10in	54101	
	8x10in	58101	
(Thin)	2x2in	52021T	
	4x5in	54051T	
	4x10in	54101T	
IFWM	3cc	FWD301	
PMX	2x2cm	-	607-005-220
	3x3cm	-	607-005-330
	4x4cm	607-001-440	607-005-440
	5x5cm	-	607-005-550
	6x6cm	607-001-660	607-005-660
	8x8cm	607-001-880	607-005-880
	8x12cm	607-001-812	607-005-812
	10x12cm	607-001-112	607-005-125
	10x25cm	607-001-125	-
	20x25cm	607-001-225	607-005-225

